



 **MEDIREX**

Positive Patient ID  
& Process Solutions

Laminated ID Bands for  
Laser Printers

# The Laminated Laser Wristband

## BARCODE READY

The Laminated Laser Wristbands are ideal for environments that require a laminated seal to preserve patient information against fading and smudging. The self-laminated design with tamper-evident adhesive closure allows for easy application. These wristbands are compatible with laser and ink jet printers, and suitable for all hospital visits up-to 7 days. The laminated laser wristbands are available in adult and infant wristband sizes, a variety of laser label quantities on the same sheet, and various colour selections.

### Patient Safety

- Laminated seal to resist moisture, alcohol, and hand sanitizer
- Hypoallergenic, latex-free, phthalate-free, and single-use disposable

### Durability

- Developed with a tamper-evident adhesive closure
- Recommended for visits lasting up-to 7 days

### Printer Compatibility


- Suited for all laser printer manufacturers
- Compatible with ink jet printers

### Easy Application

- Wristband application instructions highlighting a 4-step self-laminating process for easy application

### Highly Versatile

- Available in adult and infant wristband sizes, a variety of laser label quantities on the same sheet, and colours suited for your unique needs



# Laminated Laser Adult Wristband with 8-Labels

**BARCODE READY**

The Laminated Laser Adult Wristband with 8-Labels is the ideal solution for low acuity visits or short stays by reducing the wastage of labels when a high label count is not required. The laminated seal preserves patient information against fading and smudging, while the tamper-evident adhesive closure is tested and proven to be durable for hospital visits up-to 7 days. A 4-step self-laminating process provides easy wristband application and trouble-free operator experience. These wristbands are also hypoallergenic, latex-free, phthalate-free, and single-use disposable to further enhance the patient's safety.

Item No.	Colour
7505	<input type="radio"/> Solid White
7505R	<input checked="" type="radio"/> Red Border

### Dimensions (L x H)

Overall Size: 279MM X 147MM [11" X 5.8"]  
 Wristband: 267MM X 20MM [10.50" X 0.75"]  
 Label (8): 64MM X 25MM [2.5" X 1.0"]

**Packaging String: 1CASE/4PACKAGES/250EACH**

For more personalized solutions, please contact Medirex Sales (sales@medirex.com).

**FURTADO, SARA**  
 DOB: 06/26/86 AGE: 35Y  
 Unit #: 12345678

**STEP 1** Print ID band and labels on laser printer.  
**STEP 2** (See below)  
**STEP 3** Turn sheet over and lift ID band off sheet; start from the right and pull to the left.  
**STEP 4** Fold form over imaged area to create ID band.

**SELF LAMINATING WRISTBAND**

**MEDIREX**  
 The Laminated Wristband - 7505  
[www.medirex.com](http://www.medirex.com) | [orders@medirex.com](mailto:orders@medirex.com)

**7505**

Licensed under one or more of the following U.S. Patent 6,653,472, and LaserBand U.S. Patents 6,510,634,881; 6,067,739; 6,000,160; and others pending in the U.S. and Canada

# Laminated Laser Adult Wristband with 18-Labels

**BARCODE READY**

The Laminated Laser Adult Wristband with 18-Labels is recommended for extended hospital stays and ideal for environments requiring additional labels. The laminated seal preserves patient information against fading and smudging, while the tamper-evident adhesive closure is tested and proven to be durable for hospital visits up-to 7 days. A 4-step self-laminating process provides easy wristband application and trouble-free operator experience. These wristbands are also hypoallergenic, latex-free, phthalate-free, and single-use disposable to further enhance the patient's safety.

Item No.	Colour
7501G	Green Border
7501N	Solid White
7501R	Red Border

### Dimensions (L x H)

Overall Size: 279MM X 216MM [11.0" X 8.5"]  
 Wristband: 267MM X 20MM [10.50" X 0.75"]  
 Label (2): 127MM X 25MM [5.0" X 1.0"]  
 Label (16): 64MM X 25MM [2.5" X 1.0"]

**Packaging String: 1CASE/4PACKAGES/250EACH**

For more personalized solutions, please contact Medirex Sales (sales@medirex.com).

The image displays a sheet of 18 wristband labels arranged in a 3x6 grid. Each label is for a patient named 'FURTADO, SARA' with ID# ED00000001, PNH# 123456789-AA, GENDER F, D.O.B. 06/26/86, and AGE 35. The labels include a barcode and a QR code. The labels are shown in a grid with circular punch holes along the top edge.

**STEP 1** Print ID band and labels on laser printer.  
**STEP 2** (See below).

**STEP 3** Turn sheet over and lift ID band off sheet; start from the right and pull to the left.

**STEP 4** Fold firm over imaged area to create ID band.

**SELF LAMINATING WRISTBAND** (Left and Right sides)

**PULL RIGHT TO LEFT**

**7501N**

Push lower edge of ID band down to begin separating imaged band from sheet.

Licensed under one or more of the following U.S. Patent 6,663,472, and LaserBand U.S. Patents 6,510,634; 6,961, 6,967, 7,39; 6,000,160; and others pending in the U.S. and Canada.

# Laminated Laser Family Wristbands with 1-ID Card

**BARCODE READY**

The Laminated Laser Family Wristbands with 1-ID Card is recommended for Labour & Delivery. This family set of wristbands vastly improves patient safety by having both adult and infant wristbands on the same sheet for printing. The laminated seal preserves patient information against fading and smudging, while the tamper-evident adhesive closure is tested and proven to be durable for hospital visits up-to 7 days. A 4-step self-laminating process provides easy wristband application and trouble-free operator experience. These wristbands are also hypoallergenic, latex-free, phthalate-free, and single-use disposable to further enhance the patient's safety.

**Item No.** 7501P      **Colour**  Solid White

**Dimensions (L x H)**  
 Overall Size: 279MM X 216MM [11.0" X 8.5"]  
 Adult Wristband (2): 267MM X 20MM [10.50" X 0.75"]  
 Infant Wristband (2): 137MM X 13MM [5.625" X 0.50"]  
 ID Card (1): 82MM X 51MM [3.25" X 2.0"]

**Packaging String: 1CASE/2PACKAGES/250EACH**

For more personalized solutions, please contact Medirex Sales (sales@medirex.com).

**Adult Self-Laminating Wrist Band**

**FURTADO, SARA**  
 DOB: 06/26/86    AGE: 35Y  
 Unit #: 12345678

**STEP ONE PRINT**  
Print sheet on laser or inkjet printer.

**STEP TWO SEPARATE**  
Push down at \* and separate imaged band from sheet.

**STEP THREE FOLD**  
Begin in center of wristband - fold laminate over imaged area and smooth cut to the ends.

**STEP FOUR SECURE**  
Wrap laminated wristband around wrist. Secure with adhesive tabs.

**Adult Self-Laminating Wrist Band**

**FURTADO, SARA**  
 DOB: 06/26/86    AGE: 35Y  
 Unit #: 12345678

**Baby Self-Laminating Wrist Band**

**FURTADO, ALISHA**  
 DOB: 01/31/21  
 Unit #: 87654321

**Baby Self-Laminating Ankle Band**

**FURTADO, ALISHA**  
 DOB: 01/31/21  
 Unit #: 87654321

**Self-Laminating Identification Card**

**FURTADO, ALISHA**  
 MEMBER ID: MX0000001    SEX: F    DOB: 2021/01/31    TEL: (416) 363-9313  
 ENROLLMENT ID: LD0000001    BIRTH TIME: 08:13    TEL: (416) 363-9313  
 ATTENDING PHYS: MURPHY, SHAUN.  
 FAMILY PHYS: REZNICK, M.  
 MOTHER: FURTADO, SARA  
 EMERGENCY CONTACT: FURTADO, EDWIN

**MOTHER/FATHER/BABY**  
 SELF-LAMINATING IDENTIFICATION BANDS WITH IDENTIFICATION CARD

PLS-10SATX TAMPER EVIDENT

© 2024, LaserBand LLC. Protected by one or more of the following LaserBand U.S. Patents: #5,516,834; #5,458,881; #6,067,738; #6,060,160; #5,833,888 and others pending.



## Custom Products by request

At Medirex, we always get excited about the next challenge. We look forward to putting our years of experience, industry knowledge, and creativity to good use. Medirex is proud to develop personalized patient identification solutions that suit the unique needs of your hospital environment, are best-in-class, value-driven, and focused on safety.

For personalized solutions of wristband colours or label sizes and quantities, please reach out to our sales team:

✉ [sales@medirex.com](mailto:sales@medirex.com)

☎ 416.363.9313  
1.800.387.9848

For general inquiries please contact:

✉ [info@medirex.com](mailto:info@medirex.com)

Other ways to get in touch:

🌐 [www.Medirex.com](http://www.Medirex.com)

📘 @medirexsystems

🌐 <https://ca.linkedin.com/company/medirex>

🐦 @medirexsys

