

FURTADO, SARA
ATT. PHYSICIAN
MURPHY, SEAN
FAM. PHYSICIAN
REZNICK, M.

ED00000001

 **MEDIREX**

Positive Patient ID
& Process Solutions

ThermaSoft ID Bands for
Direct Thermal Printers

The ThermaSoft Thermal Wristband

BARCODE READY

Medirex takes a patient first approach in providing positive patient identification solutions. The ThermaSoft Thermal Wristbands are tested and proven to be #1 in the patient identification industry for thermal printers. The easy-peel closure that breaks apart to reveal a tapered adhesive end allows for trouble-free application to enhance process-of-care efficiencies. Its premium material and design improves patient comfort and safety for all hospital departments ranging from short visits to extended stays.

Superior Softness

- Material with enhanced softness to prevent discomfort for all patient skin types

Patient Safety

- Anti-wicking and tamper proof coating to protect against moisture, alcohol, and antimicrobial solutions
- Hypoallergenic, latex-free, phthalate-free, and single-use disposable

Premium Durability

- More malleable, long-lasting material tested and proven to be #1 in the patient identification industry for thermal printers
- Durable for stays up-to 14 days

Printer Compatibility

- Suited for all thermal printer manufacturers

Easy Application

- Easy-peel closure that breaks apart to reveal a tapered adhesive end for trouble-free application

Highly Versatile

- Available in adult and infant wristband sizes, a diverse range of colours, and can come with a variety of thermal labels



ThermaSoft Thermal Adult Wristband

BARCODE READY

The ThermaSoft Thermal Adult Wristband is one of the most popular products for hospitals with direct thermal printers. The easy-peel closure saves time by providing an extra edge for trouble-free removal and application. The patented material is specially designed for stays up-to 14 days and enhances patient safety by resisting moisture, alcohol, and antimicrobial solutions to preserve patient data. The ThermaSoft thermal adult wristband is hypoallergenic, latex-free, phthalate-free, single-use disposable, and designed with premium softness for all patient skin types. Various label sizes and rolls are available to complement these wristbands.


Item No.	Colour
4023-BL-D	Dark Blue Stripe
4023-BL-L	Light Blue Stripe
4023-GR	Green Stripe
4023-OR	Orange Stripe
4023-PI	Pink Stripe
4023-RE	Red Stripe
4023-WH	Solid White
4023-YE	Yellow Stripe

Dimensions (L x H)

Wristband: 279MM X 25MM [11" X 1"]
 Patient Information Area: 127MM X 25MM [5" X 1"]
 Core Size: 25MM [1"]

Packaging String: 1CASE/6BOXES/6ROLLS/200EACH

For more personalized solutions, please contact Medirex Sales (sales@medirex.com).



ThermaSoft Thermal Infant Wristband

BARCODE READY

The ThermaSoft Thermal Infant Wristband is one of the most trusted direct thermal patient identification solutions for Labour & Delivery and NICU environments. Designed with premium softness, these wristbands are very comfortable for all skin types and last for stays up-to 14 days. The anti-wicking and tamper-proof coating protects against moisture, alcohol, and antimicrobial solutions to preserve patient data and enhance patient safety. The ThermaSoft thermal infant wristband is hypoallergenic, latex-free, phthalate-free, and single-use disposable to provide a more positive patient identification experience.

Item No.	Colour
4011-BL-D	 Dark Blue Stripe
4011-BL-L	 Light Blue Stripe
4011-GR	 Green Stripe
4011-OR	 Orange Stripe
4011-PI	 Pink Stripe
4011-RE	 Red Stripe
4011-WH	 White Stripe
4011-YE	 Yellow Stripe

Dimensions (L x H)

Wristband: 152MM X 20MM [6" X 0.8125"]

Patient Information Area: 45MM X 20MM [1.8" X 0.8125"]

Core Size: 25MM [1"]

Packaging String: 1BOX/6ROLLS/350EACH

Wide Version Available


Use item codes above with -W i.e. 4011-GR-W

Wristband: 152MM X 25MM [6" X 1"]

Patient Information Area: 45MM X 25MM [1.8" X 1"]

Core Size: 25MM [1"]

Packaging String: 1BOX/8ROLLS/350EACH



For more personalized solutions, please contact Medirex Sales (sales@medirex.com).

Sample of 4011

Custom Products by request

At Medirex, we always get excited about the next challenge. We look forward to putting our years of experience, industry knowledge, and creativity to good use. Medirex is proud to develop personalized patient identification solutions that suit the unique needs of your hospital environment, are best-in-class, value-driven, and focused on safety.

For personalized solutions of wristband colours or label sizes and quantities, please reach out to our sales team:

✉ sales@medirex.com

☎ 416.363.9313
1.800.387.9848

For general inquiries please contact:

✉ info@medirex.com

Other ways to get in touch:

🌐 www.Medirex.com

📘 @medirexsystems

🌐 <https://ca.linkedin.com/company/medirex>

🐦 @medirexsys

