



 **MEDIREX**

Positive Patient ID
& Process Solutions

SuperSoft ID Bands for Laser
Printers

SuperSoft Laser Adult Wristband with 4-Labels

BARCODE READY

The SuperSoft Laser Adult Wristband with 4-Labels is the ideal solution for low acuity visits requiring the wristband only and up to four labels. Having fewer labels attached with the wristband on the same sheet helps reduce wastage. The revolutionary soft material prevents cutting or digging into the patient's skin and is durable for stays up-to 14 days. These wristbands have an anti-wicking and tamper-proof coating that withstands moisture, alcohol, and antimicrobial solutions. The SuperSoft laser adult wristband with 4-labels is also hypoallergenic, latex-free, phthalate-free, and single-use disposable to further enhance the patient's safety.

- | | |
|-----------------|---|
| Item No. | Colour |
| 8704-W | <input type="radio"/> Solid White |
| 8704-R | <input checked="" type="radio"/> Red Border |

Dimensions (L x H)

Overall Size: 279MM X 108MM [11" X 4.25"]
 Wristband: 267MM X 25MM [10.5" X 1"]
 Label (4): 64MM X 25MM [2.5" X 1"]

Packaging String: 1CASE/4PACKAGES/250EACH

For more personalized solutions, please contact Medirex Sales (sales@medirex.com).



SuperSoft Laser Adult Wristband with 20-Labels

BARCODE READY

The SuperSoft Laser Adult Wristband with 20-Labels is the ideal solution for longer stays or environments requiring an extended label count. The higher label count is optimal for high volume areas and high-acuity visits such as inpatient encounters where labels with patient information is required for hospital charts, lab results, and records. The revolutionary soft material prevents cutting or digging into the patient's skin and is durable for stays up-to 14 days. These wristbands have an anti-wicking and tamper-proof coating that withstands moisture, alcohol, and antimicrobial solutions. The SuperSoft laser adult wristband with 20-labels is also hypoallergenic, latex-free, phthalate-free, and single-use disposable to further enhance the patient's safety.


Item No.	Colour
8720-W	Solid White
8720-R	Red Border

Dimensions (L x H)

Overall Size: 279MM X 216MM [11" X 8.5"]
 Wristband: 267MM X 25MM [10.5" X 1"]
 Label (20): 64MM X 25MM [2.5" X 1"]

Packaging String: 1CASE/4PACKAGES/250EACH

For more personalized solutions, please contact Medirex Sales (sales@medirex.com).



Peel the wristband off the sheet by removing the Tab ends first. See reverse for instructions on applying wristbands to patients.

SuperSoft Laser Family Wristbands with 3-ID Cards

BARCODE READY


The SuperSoft Laser Family Wristbands with 3-ID Cards are designed for Labour & Delivery and NICU to enhance the patient identification process by having both adult and infant wristbands in one set. This vastly improves the safety for both mother and baby, and provides a better patient experience through positive patient identification. The revolutionary soft material prevents cutting or digging into the patient's skin and is durable for stays up-to 14 days. These wristbands have an anti-wicking and tamper-proof coating that withstands moisture, alcohol, and antimicrobial solutions. The SuperSoft laser family wristbands with 3-ID cards are also hypoallergenic, latex-free, phthalate-free, and single-use disposable to further enhance the patient's safety.

Item No. 8703-FAM-W
Colour Solid White

Dimensions (L x H)
 Overall Size: 279MM X 216MM [11" X 8.5"]
 Adult Wristband (2): 267MM X 25MM [10.5" X 1"]
 Infant Wristband (4): 181MM X 19MM [7.125" X 0.75"]
 ID Card (3): 86MM X 38MM [3.375" X 1.5"]

Packaging String: 1CASE/4PACKAGES/250EACH

For more personalized solutions, please contact Medirex Sales (sales@medirex.com).



Custom Products by request

At Medirex, we always get excited about the next challenge. We look forward to putting our years of experience, industry knowledge, and creativity to good use. Medirex is proud to develop personalized patient identification solutions that suit the unique needs of your hospital environment, are best-in-class, value-driven, and focused on safety.

For personalized solutions of wristband colours or label sizes and quantities, please reach out to our sales team:

✉ sales@medirex.com

☎ 416.363.9313
1.800.387.9848

For general inquiries please contact:

✉ info@medirex.com

Other ways to get in touch:

🌐 www.Medirex.com

📘 @medirexsystems

🌐 <https://ca.linkedin.com/company/medirex>

🐦 @medirexsys

