



 **MEDIREX**

Positive Patient ID
& Process Solutions

Vinyl ID Bands

The Vinyl Wristband

The Vinyl Wristbands provide an easy and effective patient identification solution for a wide variety of applications such as handwritten information, insert cards, adhering printed labels, and indicating specific patient alerts through secondary banding. The security snap closure helps minimize the risk of patients intentionally removing their wristbands. The Vinyl wristbands are suited for all hospital environments including short visits and extended stays. With three layers of comfortable and durable material, these wristbands are resistant to water, alcohol, hand sanitizers, and general cleansers.

Patient Safety

- Strong and secure snap closure to prevent tampering or reuse
- Resistant to water, alcohol, hand sanitizers, and general cleansers
- Latex and phthalate free

Durability

- Created with three layers of flexible vinyl material for stays up-to 14 days

Write-On or Imprint

- Directly write-on wristband, insert patient information card, or adhere printed label to wristband to communicate patient information

Highly Versatile

- Available in a wide range of colours and can come with a variety of labels as additional purchase
- Personalized printing of alerts onto wristband available upon request



4005

The 4005 Vinyl Adult Narrow Wristband is designed with three layers of vinyl material for extra durability and resistance to water, alcohol, hand sanitizers, and general cleansers. The snap closure ensures optimal patient safety and security. Custom pre-printed alerts are available upon request.

Packaging String: 1CASE/500EACH

Item No.

Item No.	Colour
4005-BL	● Solid Blue
4005-GR	● Solid Green
4005-OR	● Solid Orange
4005-PI	● Solid Pink
4005-PU	● Solid Purple
4005-RE	● Solid Red
4005-WH	● Solid White
4005-YE	● Solid Yellow

Dimensions (L x H)

Wristband: 259MM X 19MM [10.1875" X 0.75"]



Sample of 4005-GR

4006

The 4006 Vinyl Adult Wristband is ideal for applying handwritten patient information. The white coated area accepts permanent marker, ink pen, or multiple label sizes to clearly identify suitable information during a patient's stay. These snap closure wristbands are designed with three layers of vinyl material that is resistant to water, alcohol, hand sanitizers, and general cleansers.

Packaging String: 1CASE/250EACH

Item No.

Item No.	Colour
4006-BL	● Solid Blue
4006-GR	● Solid Green
4006-OR	● Solid Orange
4006-PI	● Solid Pink
4006-PU	● Solid Purple
4006-RE	● Solid Red
4006-WH	● Solid White
4006-YE	● Solid Yellow

Dimensions (L x H)

Wristband: 300MM X 28MM [11.875" X 1.125"]
 Patient Information Area: 95MM X 25MM [3.75" X 1.0"]



Sample of 4006-OR

4007

The 4007 Vinyl Adult Insert Card Wristband has a locking snap closure that enhances patient safety by minimizing the risk of patients intentionally removing their wristbands. The small slit on the underside of the wristband creates an opening where handwritten or pre-printed cards can be inserted. The three-layered vinyl material is resistant to water, alcohol, hand sanitizers, and general cleansers.

Packaging String: 1CASE/250EACH

Item No.

Item No.	Colour
4007-BL	● Solid Blue
4007-CL	● Solid Clear
4007-GR	● Solid Green
4007-OR	● Solid Orange
4007-PI	● Solid Pink
4007-PU	● Solid Purple
4007-RE	● Solid Red
4007-WH	● Solid White
4007-YE	● Solid Yellow

Dimensions (L x H)

Wristband: 300MM X 28MM [11.875" X 1.125"]
 Patient Information Area: 89MM X 25MM [3.5" X 1.0"]



Sample of 4007-RE

For more personalized solutions, please contact Medirex Sales (sales@medirex.com).

4004

The 4004 Vinyl Infant Narrow Wristband is designed with soft vinyl material and a re-adjustable closure. These wristbands are recommended for Labour & Delivery and NICU. Custom designs and colours are available upon request.


Packaging String: 1CASE/250EACH

Item No. Colour

- 4004-WH Solid White
- 4004-WH-TB Solid White with Teddy Bear Prints

Dimensions (L x H)

Wristband: 215MM X 13MM [8.5" X 0.5"]
 Patient Information Area: 54MM X 13MM [2.125" X 0.5"]



Sample of 4004-WH

4008

The 4008 Vinyl Adult and Infant Insert Card Wristbands come in a 1-part mother and 2-part infant wristbands set. The small slit on the underside of the wristband creates an opening where handwritten or pre-printed cards can be inserted. Designed with soft vinyl material and a re-adjustable closure, these wristbands are recommended for Labour & Delivery and NICU. Custom designs, different part combinations, and colours are available upon request.


Packaging String: 1CASE/100SETS

Item No. Colour

- 4008-WH Solid White
- 4008-WH-TB Solid White with Teddy Bear Prints

Dimensions (L x H)

Adult Wristband: 280MM X 28MM [11.0" X 1.0"]
 Infant Wristband (2): 185MM X 13MM [7.5" X 0.5"]
 Patient Information Area (Adult): 57MM X 25MM [2.25" X 1.0"]
 Patient Information Area (Infant): 32MM X 13MM [1.25" X 0.5"]



Sample of 4008-WH

4016

The 4016 Vinyl Infant Insert Card Wristband has a locking snap closure that enhances patient safety and security. The small slit on the underside of the wristband creates an opening where handwritten or pre-printed cards can be inserted. The three-layered vinyl material is resistant to water, alcohol, hand sanitizers, and general cleansers.


Packaging String: 1CASE/500EACH

Item No. Colour

- 4016-BL Solid Blue
- 4016-GR Solid Green
- 4016-OR Solid Orange
- 4016-PI Solid Pink
- 4016-PU Solid Purple
- 4016-RE Solid Red
- 4016-WH Solid White
- 4016-YE Solid Yellow

Dimensions (L x H)

Wristband: 170MM X 22MM [6.625" X 0.875"]
 Patient Information Area: 51MM X 16MM [2.0" X 0.625"]



Sample of 4016-BL

For more personalized solutions, please contact Medirex Sales (sales@medirex.com).

Custom Products by request

At Medirex, we always get excited about the next challenge. We look forward to putting our years of experience, industry knowledge, and creativity to good use. Medirex is proud to develop personalized patient identification solutions that suit the unique needs of your hospital environment, are best-in-class, value-driven, and focused on safety.

For personalized solutions of wristband colours or label sizes and quantities, please reach out to our sales team:

✉ sales@medirex.com

☎ 416.363.9313
1.800.387.9848

For general inquiries please contact:

✉ info@medirex.com

Other ways to get in touch:

🌐 www.Medirex.com

📘 @medirexsystems

🌐 <https://ca.linkedin.com/company/medirex>

🐦 @medirexsys

